Introduction

01

This project focuses on creating a higher quality method of fingerprinting that reveals more visible minutiae.

Problem

Issues with current fingerprinting methods

Traditional florescent dyes and powders:

expensive, difficult to synthesize, and aren't environmentally friendly

Small range of florescence so multiple powders are needed for different surfaces increasing variability and room for error.

Currently 22 known instances of incorrect court rulings based on fingerprint evidence and studies suggest there could be thousands more (abcnews.go.com).

Background

- <u>Carbon dots</u> are a carbonbased nanomaterial which are fluorescent and have a low toxicity, making them a good candidate for fingerprinting
- Cyanoacrylate fuming is one of many used methods of fingerprint analysis which involves using a cyanoacrylate glue to fume fingerprints
 - 1 step fluorescent dye/powder applied before fuming
 - 2 step fluorescent dye/powder applied after fuming
- Minutiae characteristics of fingerprints used to compare and classify fingerprints




Image of carbon dot florescence spectrum taken from advancedsciencenews.com




Figure of fingerprint minutiae taken from (Azman, Mahat, Wahab et al., 2019)

02_


Goals


Research Question

How can using carbon dots in conjunction with cyanoacrylate fuming create higher quality latent fingerprints?


Engineering Goals

- OT Create prints that fluoresce under 625nm, 530nm, and 465nm
- O2 Create prints with at least 5 more minutiae on average than just cyanoacrylate prints
- O3 Create prints with at least 5 more minutiae on average than just carbon dot dusted prints









Conclusion

P-Value = 0.0026 Between Control Group 1 and Experimental Group. 52 more minutiae on average.

Discussion

- · Between control Group 2 & the Experimental Group the novel procedure still produced more than 13 visible minutiae on average exceeding my engineering goal of 5 or more visible minutiae
- 39.5 more visible minutiae on average between Control Group 1 vs 2
- The NFIQ score for my novel prints was higher on average than both control groups.
 - o Higher under red and blue light when compared to Control Group 2



under FpMV software

Limitations

- · Camera quality & angle of photo
- Pressure of fingerprint deposit
- Substances available to test (couldn't test tunability) & amount of carbon dots limited
 - o limited materials carbon dots were tested on

10

Future Application

01

Test under a wider range of wavelengths using a florescent microscope

Test between wavelengths of 450-650 nm

Apply on a larger scale among forensic scientists

· Reach out to Washington State crime lab

Test florescence tunability

Coat fingerprints in different substances and test for patterns in quantum yield of florescence

Run through AFIS program

• Automated Fingerprint Identification System




Image of florescent quenching of carbon dots by a range of compounds taken from Verhagen and Kelarakis 2020.




Image of carbon dot capabilities created based on graphics from (Liu and Yang 2020) and (Sabashini, Panja, Nandini 2021)

leferences



Account Co. Co. (Co.) Convenient immercration institute in the immercratic Account (a present account for the immercratic Account for the immercratic Account (a present account for the immercratic Account for the immercratic Account for the immercratic Account (a present account for the immercratic Account Account

The control of the co